



Lab furniture



INSTRUMENT TABLES HPLC SERIES

FLEXIBLE - MOBILE - ROBUST

The HPLC table series

Our laboratory furniture is characterised by **robustness** and **functionality**. They offer everything needed in the modern analytical laboratory.

With the tendency towards more **flexibility** and **mobility** in the laboratory, the HPLC table series offers the perfect basis for your laboratory instruments. The design is trimmed for maximum **stability** and ensures that your laboratory instruments are optimally protected at all times.



Widths and depths

Widths

- 60cm (in stock)
- 80cm (on request)
- 100cm (on request)
- 120cm (in stock)

Depths

- 70cm (on request)
- 80cm (in stock)
- 90cm (on request)
- 100cm (on request)



Benefits

Flexible

- Many different widths, depths and heights available
- Electrically height adjustable
- Modular design

Mobile

- Ball bearing swivel castors with brakes
- Secure stand
- Very good running characteristics

Robust

- 250kg load capacity with even load distribution
- 150kg load capacity with uneven load distribution

Integrated cable duct

- Spacious design
- Velcro cable ties included

Integrated power distribution

- Socket strip integrated in the cable duct
- Multi-socket strips with different socket types available

Feedthroughs for cables and hoses

Fully compliant with laboratory standard 61010

Tested by an independent test laboratory

Heights

Heights for fixed height models

- 60cm (in stock)
- 70cm (on request)
- 80cm (in stock)
- 90cm (on request)
- 100cm (on request)

Height ranges for height-adjustable models

- 53-83cm -> 30cm travel (in stock)
- 58-93cm -> 35cm travel (on request)
- 63-103cm -> 40cm travel (in stock)
- 68-113cm -> 45cm travel (on request)
- 73-123cm -> 50cm travel (on request)



Materials

Tabletop

- TRESPA® TOPLAB® PLUS
- Solvent resistant
- Easy to clean
- Excellent resistance to aggressive chemicals and cleaning agents
- Does not promote bacterial growth

Substructure and chassis

- Mechanically robust steel construction with powder-coated surface

Table feet

- Aluminum



Available accessories

Drawer modules



Drawer module with soft-stop pull-outs.

- Front panel height: 119mm
- Drawer height: 68mm

Waste bottle tray



Tray with removable anti-slip mat, e.g. for waste bottles.

- Depth: 23cm
- Width: Table width - 20cm
- Tub height: 5cm

Fixed floor plate



Shelf for e.g. a PC or accessories.

- Material: TRESPA® TOPLAB® PLUS



IMPROVE YOUR SCIENCE



SONATION GMBH
ALTE SCHULSTRASSE 39
DE-88400 BIBERACH

TEL.: +49 (0) 7351 80093-0
E-MAIL: MAIL@SONATION.COM
WWW.SONATION.COM