

Sonation GmbH

for a silent lab.

Oil change

quick - easy - clean



www.sonation.com

Sonation GmbH

for a silent lab.



Oil change

quick - easy - clean

The **Sonation oil change kit** exchanges the oil drain plug of the pump by a **comfortable adapter with valve**. By opening the valve the complete oil can be drained into the supplied **oil drip sump**.

Therefore every oil change is a clean and easy process.

Oil change kit (OC-Std)

- 1 comfortable adapter with valve
- 1 white plastic hose
- 1 oil drip sump

Oil change kit for the OrbiTrap® (OC-Orbi)

- 2 comfortable adapters with valve
- 2 white plastic hoses
- 1 oil drip sump

www.sonation.com